

anti-TNF-alpha (Infliximab)

Cat# A510019

Upon receipt, store at -80°C. Avoid repeated freeze.

INFORMATION

Product Name	anti-TNF-alpha (Infliximab)
Cat. No.	A510019
Size	1mg/5mg/20mg
Isotype	Chimeric IgG1
Description	Infliximab biosimilar is a chimeric IgG1k monoclonal antibody (composed of human constant and murine variable regions) specific for human tumor necrosis factor-alpha (TNF α).
Species Reactivity	Human
Cas No	170277-31-3
Host	CHO cells
Clonality	Monoclonal
Conjugate	Non-conjugated
Tested applications	Functional Assay
Molecular Weight	149 kDa
Purity	>95%; < 0.75 EU/mg as determined by the LAL method
Storage buffer	0.01 M phosphate buffered saline (PBS) pH 7.2
Storage instruction	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at -80°C. Avoid repeated freeze/thaw cycles.
Alias	Remicade, Remsima, Inflectra

PRODUCT USE LIMITATION

These products are intended for research use only.