

Human IgM-Myeloma, Purified

Cat# P0101

Store at 2-8°C.

INFORMATION

Product Name	Human IgM-Myeloma, Purified	
Product Type	Native Protein	
Cat NO.	P0101	
Size	1mg	
Figure	 <p style="margin-left: 100px;">Lane 1. Human IgM-Myeloma, Purified Lane 2. Molecular standard</p>	
Applications:	Western Blot and ELISA	
Storage	2-8°C	
Storage Buffer	50mM TRIS, pH 8.0, 0.2 M NaCl, 0.05% (w/v) Sodium Azide	
Purity	>95% (SDS-PAGE)	
Source	Human Serum/Plasma was obtained from myeloma donors of US origin	

PRODUCT USE LIMITATION

These products are intended for research use only.