



# IL-8 mouse Monoclonal Antibody (13F8)

Cat# A22506PI

Upon receipt, store at -20°C . Avoid repeated freeze.

## INFORMATION

<b>Product Name</b>	IL-8 mouse Monoclonal Antibody (13F8)
<b>Cat. No.</b>	A22506PI
<b>Size</b>	50 µl ,100 µl
<b>Product type</b>	Primary Antibody
<b>Species Reactivity</b>	Hu, Ms, Rb
<b>Immunogen</b>	Synthetic C-TERM Peptide of IL-8
<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Tested applications</b>	IHC-p: 1:100-1:200.
<b>Conjugation</b>	Non-Conjugated
<b>Purification</b>	The antibody was affinity-purified from mouse ascites by affinity chromatography using specific immunogen.
<b>Concentration</b>	1 mg/mL
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Image</b>	 <p>Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using IL-8 Mouse mAb diluted at 1:200.</p>



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using IL-8 Mouse mAb diluted at 1:200.

=

### **PRODUCT USE LIMITATION**

These products are intended for research use only.