

20X PBS Buffer, sterile

Cat# C1037

Store at room temperature, maintain sterility

INTRODUCTION

PBS (phosphate buffered saline) is the most widely used balanced salt solution in biological research, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting, and preparing reagents. It is diluted to a 1X working concentration

Product information

Form	Clear, colorless liquid
Product size	1L
pH range	7.3-7.5
Inorganic Salts	No Calcium, No Magnesium
Phenol Red Indicator	No Phenol Red

Product components

Diluted to a 1X working concentration, contains	
NaCl	137mM
KCl	2.7mM
Na ₂ HPO ₄	8mM
KH ₂ PO ₄	2mM

PRODUCT USE LIMITATION

These products are intended for research use only.