

anti-VEGF (Bevacizumab)

Cat# A510013

Upon receipt, store at -80°C. Avoid repeated freeze.

INFORMATION

Product Name	anti-VEGF (Bevacizumab)
Cat. No.	A510013
Size	1mg/5mg/20mg
Isotype	Human IgG1
Description	Bevacizumab biosimilar is vascular endothelial growth factor directed antibody. Bevacizumab biosimilar is a recombinant humanized monoclonal IgG1 antibody that contains human framework regions and murine complementarity-determining regions. Bevacizumab biosimilar has an approximate molecular weight of 149 kDa.
Species Reactivity	Human
Cas No	216974-75-3
Host	CHO cells
Clonality	Monoclonal
Conjugate	Non-conjugated
Tested applications	Functional Assay
Molecular Weight	149 kDa
Purity	>95%; < 0.75 EU/mg as determined by the LAL method
Storage buffer	0.01 M phosphate buffered saline (PBS) pH 7.2
Storage instruction	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at -80°C. Avoid repeated freeze/thaw cycles.
Alias	rhuMab VEGF, Avastin

PRODUCT USE LIMITATION

These products are intended for research use only.