


β-actin Monoclonal Antibody(5B7)

Cat# A20250PI

Upon receipt, store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	β-actin Monoclonal Antibody(5B7)	
Cat. No.	A20250PI	
Size	100μg/50μg	
Product type	Primary antibody	
Species Reactivity	Human,Rat,Mouse,Mk,Dg,Ch,Hamster,Rb,Insect	
Host	Mouse	
Concentration	1 mg/ml	
Clonality	Monoclonal	
Tested applications	WB,IHC-p,IF	
Application	WB:1/2000-1/5000;IHC:1/200;Optimal dilutions should be determined by the end user	
Molecular Weight	43KD	
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.	
Storage buffer	PBS,pH 7.4,containing 0.02% sodium azide as Preservative and 50% Glycerol	
Storage instruction	-20°C.Do not aliquot the antibody,For research use only	
Image		Western blot analysis of 1 HEK293 2 SW480 3 HEPG2 4 MCF-7 5 mouse brain 6 Rat brain 7 HeLa 8 A549 lysates, primary antibody was diluted at 1:5000, 4° over night, secondary antibody(cat : RS23910) was diluted at 1:10000, 37° 1hour.

PRODUCT USE LIMITATION

These products are intended for research use only.