

HSP70 Polyclonal Antibody




Cat#: A20287PI

Store at -20°C. Avoid repeated freeze.

INFORMATION

Product Name	HSP70 Polyclonal Antibody
Cat. No.	A20287PI
Size	100µg/50µg/20µg
Reactivity	Human
Gene name	HSPA1A
Protein name	Heat shock 70 kDa protein 1A/1B
Human gene Id	3303/3304
Human swissprot	P08107
Immunogen	Synthesized peptide derived from human HSP70 around the non-acetylation site of K246.
Specificity	HSP70 Polyclonal Antibody detects endogenous levels of HSP70 protein only when acetylation at K246.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Rabbit
Applications	IF: 1:50-200 Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Concentration	1 mg/ml
Molecular Weight	70KD

PRODUCT IMAGES



PRODUCT USE LIMITATION

These products are intended for research use only.